



# T02 GD-OES depth profiling

## Principle

- Ar sputtering of negative electrode surface to current collector
- Detection of sputtered elements by optical emission spectroscopy (incl. Li, Si)
- Measurements at different parts of a sample due to small measuring spot ( $\varnothing$  2.5mm)



## Example: Depth profile of a Si/C electrode from a commercial 18650 cell



## Information from this method

- Fast depth profile of elements from anode surface to current collector for elements included in the calibration
- Information on aging mechanisms from Post-Mortem analysis
- More possibilities currently under development

## References for this method

- 1) SEI growth: T. Waldmann et al., Journal Electrochem. Soc. 162 (2015) A1500.
- 2) Li deposition: N. Ghanbari et al., J. Phys. Chem. C 120 (2016) 22225.
- 3) Si content of negative electrodes: K. Richter et al., Journal Electrochem. Soc. 165 (2018) A3602.