

# T02 GD-OES depth profiling - ZSW



**TEESMAT**

Open Innovation Test Bed for Electrochemical  
Energy Storage Materials

## Principle

- Ar sputtering of negative electrode surface to current collector
- Detection of sputtered elements by optical emission spectroscopy (incl. Li, Si, C, O, P, Cu, ...)
- Measurements at different parts of a sample due to small measuring spot ( $\varnothing$  2.5mm)

## Why is it useful ?

- Semi-quantitative detection of SEI growth/Li plating on graphite anodes
- Detection of Cu
- Determination of Si content in anodes


## How it works

Destructive analysis of anodes taken from Post-Mortem analysis of production line

## What can be seen

Depth profile of elemental distribution for anodes from electrode surface to current collector

**What kind of sample ?** Anodes,  $\sim 1\text{cm}^2$ , reproduction required



**Investigation time-scale** : days

**Maturity level** : advanced



Grant Agreement  
No 814106



# Example 1: SEI growth on graphite anodes




TEESMAT

Open Innovation Test Bed for Electrochemical Energy Storage Materials



GD-OES depth profile of Graphite anode after cycling at 45°C



After formation:

- Very thin initial SEI (<0.5μm) mostly on anode surface as indicated by Li peak
- Li consumption = initial capacity loss during formation

After long-term cycling (see Figure):

- SEI-thickness increased (~2μm) compared to initial SEI
- SEI growth mostly on anode surface
- Li consumption in SEI explains most of the capacity loss during aging



## Example 2: SEI growth on Si/graphite anodes




TEESMAT


Open Innovation Test Bed for Electrochemical Energy Storage Materials



GD-OES depth profile of Si/Graphite anode after formation/cycling



Raman spectroscopy of Si/Graphite anode surface



- Li peak indicates SEI formation mostly on Si/graphite anodes surface
- No Si peak before formation
- reduction of peak by washing with DMC
- Raman spectroscopy indicates Li silicate formation



# Example 3: Li deposition on graphite anodes




**TEESMAT**


Open Innovation Test Bed for Electrochemical Energy Storage Materials



GD-OES depth profile of Li plating on graphite anode after 18 cycles at 0°C / 0.5C



Anode potentials vs. Li/Li<sup>+</sup> measured in pouch full cells with Li reference electrode



Thermodynamic condition for Li deposition:  
Anode potential < 0V vs. Li/Li<sup>+</sup>  
(Li<sup>+</sup> + e<sup>-</sup> → Li)

T. Waldmann et al., J. Electrochem. Soc. 163 (2016) A1232.

- Li peak shows Li deposition mostly on graphite anode surface after aging
- Consistent with simulations on electrode level

S. Hein, A. Latz, Electrochim. Acta 201 (2016) 354.

# Example 4: Cu dissolution & re-deposition




**TEESMAT**


Open Innovation Test Bed for Electrochemical Energy Storage Materials



GD-OES depth profile of Re-deposited Cu on graphite anode after over-discharge to 0V



Anode potentials vs. Li/Li<sup>+</sup> measured in pouch full cells with Li reference electrode



Thermodynamic condition for Cu dissolution:  
 Anode potential > 3.56V vs. Li/Li<sup>+</sup>  
 (dissolution:  $\text{Cu} \rightarrow \text{Cu}^+ + \text{e}^-$   
 re-deposition:  $\text{Cu}^+ + \text{e}^- \rightarrow \text{Cu}$ )

- Li peak shows Cu re-deposition on graphite anode surface after over-discharge to 0V
- No Cu dissolution for cell opened at 2V
- According to Gibbs Phase Rule:  
 Constant anode potential at 3.56V vs. Li/Li<sup>+</sup>  
 → dissolution as Cu<sup>+</sup> (Cu<sup>2+</sup>: 3.38V vs. Li/Li<sup>+</sup>)